

## Bunnell Life Pulse High Frequency Ventilator Circuit

### WARNING:

- The integrated pneumatic backup (250 mL/min.) should not be used with the Bunnell Life Pulse as ventilator flow rates are normally below the recommended ventilator flows.
- Place the Bunnell Life Pulse in standby prior to suctioning the patient to avoid NO delivery transiently exceeding the set dose by up to 30 ppm. Press ENTER to reestablish ventilation as soon as the catheter is removed from the airway. This will limit the extent of over delivery above the NO set dose.

### Caution:

- If set dose is below five ppm and the Servo pressure is 2.0 psig. or less, this will result in flow rates outside of the specification of the injector module and fluctuating NO values may result.
- A one-way valve should be placed between the injector module and the humidifier chamber to prevent water from backing up into the injector module if the Life Pulse is either put into Standby or cycled OFF.
- There are higher pressures in the breathing circuit than normal; use only parts provided in disposable package #50248 and tightly secure all connections.





Figure 3-24 Example: Bunnell Life Pulse Ventilator Diagram

### Connection Instructions:

1. Connect the sample Tee as shown in [Figure 3-25](#).
2. Connect the injector module as shown in [Figure 3-26](#). The one-way valve prevents water from backing up into the injector module if the Life Pulse is either put into standby or cycled OFF.


## Connecting INOmax DS<sub>IR</sub> Sample Tee to the Bunnell Life Pulse Circuit



1. From Patient Box
2. Cut Green tube at midpoint  
(approximately six inches from the Life Port Adapter)
3. Patient Gas Sample Line with Nafion
4. Insert Sample Tee
5. Life Port Adapter
6. Endotracheal Tube

Figure 3-25

## Connecting INOmax DS<sub>IR</sub> Injector Module to the Bunnell Life Pulse Circuit



- |                                     |   |
|-------------------------------------|---|
| 1. Injector Module Electrical Cable | 6. Injector Module                      |
| 2. NO/N <sub>2</sub> Injector Tube  | 7. 15M X 4.5 mm I.D. Adapter            |
| 3. Gas Out Tube from Vent           | 8. Three cm Piece of Green Gas Out Tube |
| 4. 15M X 4.5 mm I.D. Adapter        | 9. One-Way Valve                        |
| 5. 22M/15F X 22M/15F Adapter        | 10. Green Gas Out Tube to Humidifier    |

Figure 3-26