


THE IMPORTANCE OF SERVO PRESSURE

- Servo Pressure = driving pressure that automatically regulates flow.
- Servo Pressure changes as lung volume or mechanics change.



Servo Pressure changes: early warning of changes in patient condition.

SERVO PRESSURE INCREASES WITH:

- Improving compliance or resistance
- Leak around ETT
- Tubing leak

SERVO PRESSURE DECREASES WITH:

- Worsening compliance or resistance
- Obstructed ETT
- Tension Pneumothorax
- Right mainstem intubation
- Patient needs suctioning

Charting Servo Pressure simplifies patient management decisions.