

# LifePort™ Endotracheal Tube Adapter

## Jet Ventilation Without Reintubation



### Pressure Monitoring Tube

- Provides accurate pressure regulation and reliable pressure monitoring

### 15 mm Connection

- Standard connection to conventional ventilator

### Jet Port

- Breaths from the Life Pulse ventilator are redirected toward the lungs through a jet nozzle

### Jet Port Cap

- Seals the port when not in use

### E.T. Tube Connection

- Sizes available to fit most infant endotracheal tubes

Simply replace the 15 mm adapter on the standard ET tube with the LifePort adapter, connect the Life Pulse circuit, and start ventilating.

### Life Pulse®




High Frequency "Jet" Ventilator

### Ordering Information:

|      |  |       |
|------|--|-------|
| 937  | HFJV Patient Circuit Kit                   | 2 ea  |
|      | (Includes 2 circuits, 4 LifePort adapters) |       |
| 9025 | LifePort 2.5 mm I.D.                       | 10 ea |
| 9035 | LifePort 3.5 mm I.D.                       | 10 ea |
| 9045 | LifePort 4.5 mm I.D.                       | 10 ea |

1-800-800-HFJV  
[www.bunl.com](http://www.bunl.com)

 **bunnell**  
INSPIRED INFANT CARE  
436 Lawndale Drive  
Salt Lake City, Utah 84115