

“WhisperJet” Patient Box

Quiet, Accurate and Reliable



Sealed Enclosure

- Dramatically reduces noise and keeps fluids out

Pressure Transducer

- Positioned to minimize potential fouling

Pinch Tubing Capture

- Designed to hold the pinch tubing more securely

Tubing Connections

- Now located on the ends of the box to protect against kinking

Flexible Cable

- Reduced stiffness allows for hassle-free placement

Life Pulse



High Frequency “Jet” Ventilator

The “Whisper Jet” Patient Box makes the Life Pulse High-Frequency Ventilator, a uniquely effective mode of ventilation. It allows the inhalation valve and pressure transducer to be placed close to the patient. This close proximity makes the small high velocity breaths more effective and facilitates fast and accurate pressure monitoring. The result is a system that responds to changes in your patient, moment to moment.

1-800-800-4358

www.bunl.com


INSPIRED INFANT CARE

436 Lawndale Drive
Salt Lake City, Utah 84115