

# “WhisperJet” Patient Box

## Quiet, Accurate and Reliable



### Sealed Enclosure

- Dramatically reduces noise and keeps fluids out

### Pressure Transducer

- Positioned to minimize potential fouling

### Pinch Tubing Capture

- Designed to hold the pinch tubing more securely

### Tubing Connections

- Now located on the ends of the box to protect against kinking

### Flexible Cable

- Reduced stiffness allows for hassle-free placement

## Life Pulse



High Frequency “Jet” Ventilator

The “Whisper Jet” Patient Box makes the Life Pulse High-Frequency Ventilator, a uniquely effective mode of ventilation. It allows the inhalation valve and pressure transducer to be placed close to the patient. This close proximity makes the small high velocity breaths more effective and facilitates fast and accurate pressure monitoring. The result is a system that responds to changes in your patient, moment to moment.

1-800-800-4358

[www.bunl.com](http://www.bunl.com)

 **bunnell**  
INSPIRED INFANT CARE

436 Lawndale Drive  
Salt Lake City, Utah 84115