

Closed Suction System Options

AirLife (provided by CareFusion, 800-231-2466): this system works well with the Bunnell LifePort adapter. The AirLife system comes in various sizes:

<u>Cat. No.</u>	<u>French Size</u>
50-3505	5 Fr
50-3506	6 Fr
50-3507	8 Fr

[Click for AirLife webpage](#)

Halyard: this Closed Suction Catheter, Manifold Elbow, Neonatal/Pediatric system uses Ballard technology. The sheathed catheter catalogue number and french sizes are listed below:

<u>Cat. No.</u>	<u>French Size</u>
201	5 Fr
200	6 Fr
205	8 Fr

[Click for Halyard webpage](#)

B&B Medical Technologies (800-242-8778): the NEO₂-Safe is a circuit “Y” connector with a duckbill valve built into it, so you can suction using traditional suction catheters without losing pressure in the ventilator circuit.

[Click for B&B webpage](#)