

IMPROVED PRODUCT ANNOUNCEMENT

Re-designed “FUNNEL” LifePort™ ETT Adapters for the LifePulse™ HFJV

Dear Bunnell Customer,

March 22, 2018

Bunnell Inc. is proud to provide new LifePort™ ETT Adapters with built-in *funnels* for use with the LifePulse® High-Frequency “Jet” Ventilator. Bunnell re-engineered the Adapters after many clinicians asked for a better design to allow for an easier and smoother insertion of a suction catheter, without adding dead space.



LifePort Funnel Adapter

- Funnel LifePort Adapters are now included in the HFJV Patient Circuit Kit (Cat#938), with two adapters in each size: 2.5mm, 3.0mm, and 3.5mm.
- Funnel LifePort Adapters will be available in late May in boxes of 10 in these sizes: 2.5mm, 3.0mm, 3.5mm, and 4.5mm. Notification will be sent when they are available.

Existing LifePort Adapters can still be used. To distinguish old adapters from the “funnel” adapters, please note that “FUNNEL” is added to the LifePort packaging, found only inside the Kits (Cat#938):

